

## **SINGLE LIMB CIRCUIT WITH WHISPER VALVE: Measures/Targets**

### **Expiratory Tidal Volumes**

**MODE: PS ( pressure added to whatever the peep is set)**

**PS: (Pressure Support)** Is a pressure targeted ventilation mode delivering pressure supported spontaneous breaths. (**Safety Vt optional**) **FREQUENTLY USED MODE**

Ordered Settings without Safety Tidal Volume: RR(minimum back-up rate), PEEP, PS, Ti min, Ti max

**PS SVT:** Ordered Settings with Safety Tidal Volume: RR (minimum back-up rate), PEEP, PS(now a minimum pressure) ,I-Time min, I-Time max, Tidal Volume, PS max

- With safety Volume on the vent attempts to increase the ps until the targeted exvt is obtained. Slower in resonance to exvt than the servo's (prvc or vol support modes ).
- The set rate is a back up minimum rate as the target volume and insp time functions in the same maneuver when the rate is initiated using flow termination within the allowable times.

#### **RANGES/DEFAULT VALUES**

**RR:** 5-80 bpm [20bpm]

**PEEP:** 2 -25cmH<sub>2</sub>O [5cmH<sub>2</sub>O]

**PS (above PEEP):** 2-48cmH<sub>2</sub>O [7cmH<sub>2</sub>O]

**Cycle %:** 5-90%

**Trigger (pressure):** Very Low to Very High [Very High]

**Rise Time(msec):** 150 - 900msec [200msec] (**cannot exceed**  
 **$\frac{2}{3}$  of I-time max**)

**I-time min:** 0.2-4.0 sec [0.2sec]

**I-time max:** 0.3-4.0 sec [0.8sec]

**Safety Tidal Volume (**optional**):** 50-500ml

**PS max(mandatory with safety tidal volume):** 2-48cmH<sub>2</sub>O

---

## **CPAP:**

(**Continuous Positive Airway Pressure**) mode delivers a constant level of pressure during inspiration and expiration.

#### **RANGES/DEFAULT VALUES:**

**CPAP (PEEP):** 3-20cmH<sub>2</sub>O [5cmH<sub>2</sub>O]

**Trigger (pressure):** Very low to Very High [Very High]

## **(S)T: (Spontaneous Time)**

Is a pressure targeted, flow cycled, ventilation mode delivering pressure supported spontaneous breaths at a set back-up rate with a **(Safety Vt optional)** **FREQUENTLY USED MODE**

Ordered Settings without Safety Tidal Volume: RR(back-up), EPAP, IPAP(PIP), I-Time min, I-Time max

Ordered Settings with Safety Tidal Volume: RR, EPAP, IPAP(PIP), I-Time min, I-Time max Tidal Volume, IPAP Max (pip max as set)

### **RANGES/DEFAULT VALUES**

**RR:** 5 - 80BPM [20bpm]

**EPAP (PEEP):** 2 - 25 cmH<sub>2</sub>O [5cmH<sub>2</sub>O]

**IPAP:** 4 - 50 cmH<sub>2</sub>O [12cmH<sub>2</sub>O]

**Trigger (pressure):** Very low to Very High [Very High]

**Rise Time (msec):** 150-900 [200 msec] **(cannot exceed  $\frac{2}{3}$  of I-time max)**

**I-time min:** 0.2 - 4.0 sec [0.2sec]

**I-time max:** 0.3 - 4.0 sec [0.8sec]

**Cycle%:** 5-90% [25%]

**Safety Tidal Volume (optional):** 50-500ml[50ml]

**I-PAP max (mandatory with Safety Tidal Volume):** 4.0-50 cmH<sub>2</sub>O

---

**P(A)C:** Is a bilevel ventilation mode delivering pressure-controlled mandatory breaths **(Safety Vt optional)**

Ordered Settings without safety Tidal Volume: RR, EPAP, IPAP, I-Time

Ordered Settings with safety Tidal Volume: RR, EPAP, IPAP(minimum pip), I-Time, Safety Tidal Volume, IPAP max

### **RANGES/DEFAULT VALUES:**

**RR:** 5-80 bpm[20bpm]

**EPAP (PEEP):** 2 to 25 cmH<sub>2</sub>O [5cmH<sub>2</sub>O]

**IPAP:** 4 to 50 cmH<sub>2</sub>O [12cmH<sub>2</sub>O]

**I-time:** 0.3 to 4.0 sec[0.6sec] **(cannot exceed  $\frac{2}{3}$  of 60bpm)**

**Trigger (pressure):** Very Low to Very High [Very High]

**Rise Time (msec):** 150-900 msec[200 msec]

**Safety Tidal Volume (optional):** 50-500ml[50ml]

**I-PAP max (mandatory with Safety Tidal Volume): 4.0- 50 cmH<sub>2</sub>O**

**DISCONNECTION ALARM: (has to be set "ON" to use)**

**The disconnection alarm can be adjusted based on patient minimal leak**