

SINGLE LIMB CIRCUIT WITH ACTIVE EXPIRATORY VALVE:

The Astral Device **DOES NOT** support monitoring of Exhaled Volumes when used with a single limb circuit with an active expiratory valve. **Measures (targeted) Inspiratory Tidal Volumes.**

Volume modes

(A)CV (Assisted volume-controlled ventilation): Is a volume targeted ventilation mode delivery volume-controlled mandatory breaths

Ordered Settings: RR, PEEP, Tidal Volume, I-Time

APNEA Backup (Alarm only, ACV+Alarm, PACV+ Alarm)

Ordered Settings for apnea back up PIP/Tidal Volume, PEEP, RR, I-Time

RANGES/DEFAULT VALUES:

RR: 5-80bpm [20 bpm]

PEEP: 3-20 cmH₂O [5cmH₂O]

Vt: 50-300 cmH₂O [50-60cmH₂O]

PIF: (When volume breath option is set to PIF): 5-60 L/M
[10lpm]

I-time: (When volume breath option is set to Ti): 0.3-3.0sec
[0.6sec]

Flow shape %: 100(constant) 25,50,100 [100%]

Trigger (pressure): very low to very high [Very High]

V-SIMV (Volume Synchronised Intermittent Mandatory Ventilation): Is a mixed ventilation mode delivering volume-controlled mandatory breaths and pressured supported spontaneous breaths.

Ordered Settings: RR, PEEP, PS, Vt, I-Time

RANGES/DEFAULT VALUES:

RR: 5-80 bpm [20bpm]

PEEP: 3-20cmH₂O [5cmH₂O]

PS (above PEEP): Spontaneous Breaths: 2-50 cmH₂O [7cmH₂O]

Vt: 50-300 ml [50-60 ml]

PIF: (When Volume Breaths option is set to PIF) Mandatory
Breaths: 5-60 bpm [10bpm]

I-time: (When Volume Breaths option is set to I-time)
Breaths: 0.3 to 3sec [0.6sec]

Flow shape %: 100(constant) 25,50,100 [100%]
Cycle %: Spontaneous Breaths:5 to 90 [Auto]
Trigger (pressure) :Very Low to Very High [Very High]

Pressure Modes

PS: (**Pressure Support**) Is a pressure targeted ventilation mode delivering pressure supported spontaneous breaths. (**Safety Tidal Volume optional**) (**FREQUENTLY USED MODE**)
Ordered Settings without Safety Tidal Volume: RR(minimum), PEEP,PS, Ti min,Ti max
Ordered Settings with Safety Tidal Volume: RR, PEEP,PS(minimum) ,I-time max, I-time min, Tidal Volume(insp vol targeted), PS max

RANGES/DEFAULT VALUES:

RR: 5-80 bpm [20bpm]
PEEP: 2- 25cmH2O [5cmH2O]
PS (above PEEP): 2- 48cmH2O [7cmH2O]
Cycle %: 2 - 90%
Trigger (pressure):Low to Very High [Very High]
Rise Time(msec): 150 - 900msec [200msec] (**cannot exceed**
 $\frac{2}{3}$ of I-time max)

I-time min: 0.2-4.0 sec [0.2sec]

I-time max: 0.3-4.0 sec [0.8sec]

Safety Tidal Volume (optional):50-500ml inspvt

PS max(mandatory with safety tidal volume):2-50cmH2O

CPAP: (**Continuous Positive Airway Pressure**): mode delivers a constant level of pressure during inspiration and expiration.

APNEA Backup (Alarm only, ACV+Alarm, PACV+Alarm)

Ordered Settings for apnea back up PIP/Tidal Volume, PEEP, RR, I-Time

RANGES/DEFAULT VALUES:

CPAP (PEEP): 3-20cmH2O [5cmH2O]

Trigger (pressure): Very low to Very High [Very High]

P(A)CV (Assisted Pressure Controlled Ventilation): Is a pressure targeted ventilation mode delivering pressure-controlled mandatory breaths. (Safety Tidal Volume optional)

Ordered Settings without Safety Tidal Volume: RR, PEEP, P Control, I-Time

Ordered Settings with Safety Tidal Volume: RR, PEEP, Insp Tidal Volume, I-Time, P control max

Apnea Backup optional

RANGES/DEFAULT VALUES:

RR: 5- 80 bpm [20bpm]

PEEP: 3-20cmH₂O [5cmH₂O]

P Control: 2 - 50cmH₂O [7cmH₂O]

Ti(sec): 0.2-5.0 [0.6sec]

Trigger (pressure): low to Very High [Very High]

Rise Time (msec): 150-900msec [200 msec]

Safety Tidal Volume (optional): 50-500ml inspVT

P-Control max(mandatory with safety tidal volume) 2-50cmH₂O

P-SIMV (Pressure Synchronized Intermittent Mandatory Ventilation): Is a mixed ventilation mode delivering pressure-controlled mandatory breaths and pressure supported spontaneous breaths.

Ordered Settings: RR, PEEP, Pcontrol, PS, I-time

RANGES/DEFAULT VALUES:

RR: 5-80bpm [20bpm]

PEEP: 3-20cmH₂O [5cmH₂O]

P Control: 2-50 cmH₂O [7 cmH₂O]

PS (above PEEP): Spontaneous Breaths: 2-50 cmH₂O
[7cmH₂O]

I-time: 0.2-5sec [0.6sec]

Cycle %: Spontaneous Breaths: 5 - 90 [Auto]

Trigger (pressure): Very low to Very High [Very High]

Rise Time: 150 - 900 (msec) [200 msec]

DISCONNECTION ALARM: (has to be set "ON" to use)

The disconnection alarm can be adjusted based on patient minimal leak