

DUAL LIMB CIRCUIT: Measures Inspiratory and Expiratory Tidal Volumes. Exhalation Returns to device

Volume Modes

(A)CV (Assisted volume-control ventilation): Is a volume targeted ventilation mode delivery volume controlled mandatory breaths.

Ordered Settings: RR, PEEP, Tidal Volume, I-Time

RANGES/DEFAULT VALUES:

RR: 5-80 bpm [20bpm]

PEEP: 3-20cmH₂O [7cmH₂O]

I-Time: 0.3-3.0 sec [1.00]

Flow Shape: 25%-100% [100%]

Trigger (pressure): very low to Very High [very high]

V-SIMV (Volume Synchronized Intermittent Mandatory Ventilation): Is a mixed ventilation mode delivering volume controlled mandatory breaths and pressured supported spontaneous breaths.

Ordered Settings: RR, PEEP, PS, Tidal Volume, I-Time

RANGES/DEFAULT VALUES:

RR: 5-80 bpm [20 bpm]

Vt: 50-500ml [50ml]

PEEP: 3- 20cmH₂O [7cmH₂O]

PS (above PEEP): 2-50cmH₂O

I-Time: 0.3-3.0sec [1.00sec]

Flow shape: 25% to 100% [100%]

Trigger (pressure): very low to very high

Rise Time (msec): 150-900 msec [200msec]

Cycle %: 5 to Auto [5%]

Pressure Modes

P(A) CV (Assisted Pressure Controlled Ventilation) Is a pressure targeted ventilation mode delivering pressure-controlled mandatory breaths. (**Safety Tidal Volume optional**)

Ordered Settings: RR, PEEP, P-control, I-Time

Ordered Settings with Safety Tidal Volume: RR, PEEP, P-control, I-time, Tidal volume, P-control max

RANGES/DEFAULTS VALUES:

RR: 5 - 80bpm [20 bpm]

P-Control: 2 -50cmH₂O [7cmH₂O]

PEEP: 3-20cmH₂O [7cmH₂O]

I-Time: 0.3-3.0sec [0.6 sec]

Trigger (pressure): very low to very high

Rise Time (msec): 150-900 msec [200msec]

Safety Tidal Volume (optional**):**50-500ml[50ml]

P-Control max(mandatory with safety tidal volume)2-50cmH₂O [30cmH₂O]

P-SIMV (Pressure Synchronized Intermittent Mandatory Ventilation): Is

a mixed ventilation mode delivering pressured-controlled mandatory breaths and pressure supported spontaneous breaths.

Ordered Settings: RR, PEEP, PS, P-control, I-Time

Apnea Alarm (Alarm only, ACV + Alarm, PACV + Alarm)

RANGES/DEFAULTS VALUES:

RR: 5-80 bpm [20 bpm]

PEEP: 3-20cmH₂O [7cmH₂O]

P-Control: 2-50cmH₂O

I-Time: 0.2-5.0 sec [0.6 sec]

PS (above PEEP): 2- 50 cmH₂O

Rise Time (msec): 150-900 [200msec]

Trigger (pressure): very low to very high [very high]

Cycle %: 5 to auto [5%]

PS (Pressure Support): Is a pressure targeted ventilation mode delivering pressure supported spontaneous breaths (**Safety Tidal Volume optional**)

Ordered Settings without Safety Tidal Volume: RR, PEEP, PS, I-time max, I-time min

Ordered Settings with Safety Tidal Volume: RR, PEEP, PS, I-time max, I-time min, Tidal volume, PS max

RANGES/DEFAULT VALUES:

RR: 5-80 bpm [20 bpm]

PEEP: 3-20cmH₂O [7cmH₂O]

PS (above PEEP): 2-50 cmH₂O

I-time min: 0.2-4.0 sec [0.2sec]

I-time max: 0.3-4.0 sec [0.8sec]

Rise Time (msec): 150 - 900 msec [200msec] (**cannot exceed**
 $\frac{2}{3}$ of I-time max)

Trigger (pressure): very low to very high [very high]

Cycle %: 5 to auto [5%]

Safety Tidal Volume (optional): 50-500ml [50ml]

PS max (mandatory with Safety Vt): 2-30cmH₂O [30cmH₂O]

CPAP (Continuous Positive Airway Pressure): mode delivers a constant level of pressure during inspiration and expiration.

RANGES/DEFAULT VALUES:

CPAP (PEEP) 3-20cmH₂O [7cmH₂O]

Trigger (pressure): Very low to Very High [Very High]